LEICA DMIL MICROSCOPE INSTRUCTIONS MANUAL

Available languages

EN

DE

FR

Quick Links

[The Microscope and Its Components](#)

[Troubleshooting, Lamp/Fuse Replacement](#)

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Power unit

Fluorescence filter block slide

Incident-light fluorescence facility

Rear of the microscope

Ambient conditions

Installation site

Unpacking

Attaching the condensers

Stabiliser plate

Microscope with transmitted-light illumination column

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Attaching the stage expansion

Inserting the specimen brackets

Specimen stage with expansions attached

Selecting the mains voltage and connecting microscope to mains supply

Positioning the filter block slide

Assembling the fluorescence filter blocks

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Rear panel of 12 V/100 W Leica

12 V/100 W lamp holder

Lamp holders for gas discharge lamps

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Microphotography

TV adaptation

Adaptation of imaging systems on binocular photo tubes

Positioning the drawing facility

Positioning the multi-discussion facility

Microscope with drawing facility

Positioning the magnification changer

Positioning the ergo module

Switching on the 6 V/35 W halogen lamp

Adjustment specimen

Operation

Basic settings for transmitted light

Focusing on the object

DM ILB binocular tube

Adjusting tubes and eyepieces

Correction for defective vision

DM ILT trinocular tube

DM ILB tube

DM ILT tube

Locking of objectives

CORR objectives

Immersion objectives

Operation of objectives

Bright-field illumination

Adjusting the condenser

Adjusting the aperture diaphragm

Operation of transmitted light

Possible errors

Transmitted-light illumination unit with condenser

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Optimising the IMC

IMC-Components

Operation of incident light fluorescence

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Applications of condensers S 90 and S 55

Performance data of condensers

Incident-light fluorescence illumination

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Specifications of the power unit

Tools/Adjustment keys

Screw covers for unused nosepiece positions

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
[Déclaration de conformité UE](#)

[Manual Versions in Other Languages](#)

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
Leica

Copyrights

All rights to this documentation are owned by Leica Microsystems Wetzlar GmbH. Copying of text and illustrations – in full or in part – by printing, photostat, microfilm or other techniques, including electronic systems, is only permitted subject to the express written consent of Leica Microsystems Wetzlar GmbH.

The information contained in the following documentation represents the latest stage of technology and knowledge. We have composed the texts and illustrations with great care. However, as it is impossible to eliminate the risk of error completely, we cannot accept any kind of liability for the correctness of the contents of this manual. Nevertheless, we are always grateful to be notified of any errors.

The information in this manual may be altered without prior notice.

Contents

| | | | |
|--|-----------|--|-----------|
| Important notes | 6 | Technical description | 67 |
| General safety information | 7 | Objectives | 67 |
| Application | 9 | Eyepieces | 70 |
| The microscope and its components | 10 | Filters | 72 |
| Installation site | 15 | Tubes | 73 |
| Unpacking | 15 | Condensers | 75 |
| Setting up | 17 | Lamps and lamp housings | 76 |
| Setting up the options | 26 | Specifications | 78 |
| Operation | 39 | Main wear and replacement parts | 80 |
| Basic settings for transmitted light | 39 | EU Declaration of Conformity | 81 |
| Operation of objectives | 42 | | |
| Operation of transmitted light | 43 | | |
| Operation of phase contrast | 45 | | |
| Operation of Integrated Modulation Contrast (IMC) | 47 | | |
| Operation of incident light fluorescence | 51 | | |
| Care and maintenance | 56 | | |
| Troubleshooting (lamp/fuse replacement) | 58 | | |
| Storage | 66 | | |
| Packaging and transport | 66 | | |

Important notes

This manual is an important part of the Leica DMIL microscope and should be carefully read before switching on and using the instrument.

This manual contains important instructions and information with regard to operating safety and maintenance of the system. For this reason, it should be carefully kept at hand.

The manual is multi-lingual. The spiral binding makes it possible to turn the instructions in the desired language to the front.

Text symbols and their meaning:

(1.2)

Numbers in brackets, e.g. (1.2), refer to illustrations, in the example to Fig. 1, pos. 2.

→ p. 20

Numbers with an arrow, e.g. → p. 20, refer to a particular page of this manual.

Special safety information is marked with a triangular symbol in the margin and printed on a grey background.

Attention! Operating errors can damage the microscope or its accessories.

Warning of hot surface.

Explanatory note.

*

Not part of all configurations.

General safety information

This instrument of Safety Category I has been built and tested according to EN 61010-1/ IEC 1010-1, Safety Standards for Electronic Measuring Instruments, Electronic Regulators and Electronic Laboratory Instruments.

Attention!

To maintain this condition and to ensure safe operation, the user must note and adhere to the directions and warnings contained in this manual.

The mains plug may only be inserted into a grounded socket.

The protection must not be jeopardised by using an extension cable without ground conductor. Any break in the ground conductor within or outside the instrument or loosening of the ground connection can render the unit dangerous. Intentional severance is prohibited!

Attention!

Accessory devices connected to the microscope which have a separate and/or different power supply should be brought to the same ground potential by connecting them to the same grounding system. Contact the service department in case of an ungrounded mains supply .

Only fuses of the specified type and rating may be used as replacements. Never use repaired fuses or short-circuited fuse holders.

Attention!

The instruments and accessories described in this manual have been safety-tested and checked for possible hazards.

Before modifying the instrument in any way or combining it with non-Leica products not covered by this manual, please always contact your local Leica representative or the main factory in Wetzlar!

Any unauthorised interference with the instrument or use of the instrument for applications for which it is not designed will automatically void any warranty claim!

Attention!

The electrical accessories for the microscope are not protected against the penetration of water which can result in electrical shock.

Do not set up the microscope and its accessories in the immediate vicinity of a water outlet or in other places where water may penetrate into the instrument.

Attention!

Protect the microscope against excessive temperature fluctuations which may lead to condensation and thus damage the electrical and optical components.

Attention!

When using immersion oil, take care to avoid skin contact! Ask the supplier for a safety data sheet!

Attention!

Before replacing fuses or lamps, always switch off the power switch and disconnect the power cord.

Application

The Leica DM IL microscope is designed for routine examinations of cell and tissue cultures, liquids, and sediments.

Two basic microscope stands are available for biologic applications. One is the standard bright field version using the contrasting methods bright field (BF), phase contrast (Phaco), and Integrated Modulation Contrast (IMC), the other is the fluorescence microscope which also offers incident-light fluorescence in addition to the three transmitted-light contrasting techniques.

All microscope techniques and the required accessories for the Leica DM IL are described and explained in detail in the Operation section of this manual, including the associated functions and their use.

The microscope and its components

Main assemblies

The following overall views show and designate the main assemblies of the microscope and its accessories.

Fig. 1 Right side of microscope stand


Fig. 2 Left side of microscope stand



Fig. 1–3

1 Binocular photo tube DM ILT, 2 Eyepiece tube, 3 Eyepieces, 4 Tube holder, 5 Empty slide or IMC modulator, 6 Integrated lamp housing with halogen lamp 6 V/35 W, 7 Holder for filter Ø 32 mm, 8 Empty slide or modulation or phase contrast slide, 9 Aperture diaphragm, 10 Condenser S 55, 11 Transmitted-light illumination carrier, 12 Notched lever for condenser level adjustment, 13 Transmitted-light illumination column, 14 Specimen stage, 15 Objective nosepiece and objectives, 16 Brightness control, 17 Coarse control adjustment and micrometer adjustment, 18 Power switch, 19 Fluorescence filter blocks, 20 Fluorescence lamp housing, 21 Light stop, 22 BG9 Filter, 23 Stabiliser plate, 24 c-mount video adapter, 25 Tube adapter IL/L, 26 DML tubes, 27 Multi-discussion facility, 28 Condenser S 90

Fig. 3 Microscope stand with DM L tubes and discussion facility

Microscope stand

Owing to its low centre of gravity, the Leica DM IL microscope stand is very stable. When using the multi-discussion facility* or for long-term exposures in micro-photography applications, a stabiliser plate* is available to improve stability. For incident-light fluorescence applications, an incident-light axis is integrated in a second stand version.

Tube holder

The tube holder is the interface between the microscope stand and the tube. The tube holder permits the use of DM ILB and DM ILT tubes as well as the IL/L tube adapter which in turn permits the use of DM L tubes. The ergo module and the drawing module can also be directly mounted on the tube holder if the DM ILB or DM ILT tubes are used (see also Tube adapter).

Tube

The tube contains a 1x tube lense which creates the primary image in conjunction with the objective.

The binocular tube consists of a body, the binocular part, and the tube change ring.

The trinocular tube offers an additional documentation output to accommodate photo or video equipment. A switchable mirror directs the light 100% either to the eyepieces or the photo output.

DM IL tubes can be replaced and rotated.

IL/L tube adapter

The tube adapter serves to accommodate tubes from the DM L range as well as the drawing facility*, the multi-discussion facility*, the magnification changer*, and the ergo module*.

Eyepieces

The eyepiece creates a magnified, virtual image of the real intermediate image which is projected by the objective. In this process, the eyepieces works as a magnifier.

Transmitted-light illumination unit

The transmitted-light illumination facility consists of the transmitted-light illumination carrier and the transmitted-light illumination column. The transmitted-light illumination carrier contains a pre-centred, bright 6 V/35 W halogen lamp, as well as a holder for a diaphragm slide, a holder for a light filter, a condenser and an aperture diaphragm.

Lamp housing

The Leica DM IL microscope uses an integrated lamp housing with a 6 V/35 W halogen lamp.

Filters

The filters are usually employed to improve the contrast of the specimen. They are firmly mounted in a holder (Ø 32 mm).

Different filters can be inserted into the filter holder of the transmitted-light illumination unit.

Aperture diaphragm

The aperture diaphragm determines the resolution, the depth of field and the contrast of the microscope image. The best resolution is obtained when the apertures of the objective and the condenser are roughly the same.

Attention:

The aperture diaphragm in the **illumination light path** is **not** for setting the image brightness. Only the rotary brightness adjustment knob or the neutral density filter should be used for this purpose.

Condenser

The condenser is a lens system which collects the light and directs it to the specimen from the top. The condenser permits to utilise the numerical aperture in the objective.

Notched lever for condenser level adjustment

The notched lever permits the adjustment of the condenser height by moving the transmitted-light carrier. The markings on the transmitted-light column indicated the proper height to be set for the condenser used.

Specimen stages and accessories

The specimen stage is used to support the specimens to be subjected to microscopic examination. Several options are available for different objects, such as specimen controllers, stage expansions, scanning table, clamps, heating table etc.

Objective revolver and objectives

The objective revolver is designed to hold the objectives. Especially the long-range L-type objectives are able to compensate for the different thicknesses of vessel bottoms.

All microscope objectives with magnifications 2.5 to 100 can be used. All objectives from the DML and DMR range with 25 mm thread are compatible. A choice of objectives is contained in chapter "Specifications; Performance Data" or in the current objective lists that are available from your local Leica representative.

Brightness adjustment

The microscope stand contains a 6 V/35 W transformer for continuous brightness adjustment using a rotary brightness switch.

Coarse and fine focus adjustment

Coarse and fine focus adjustment permits a fast and precise adjustment of the microscopic image. Focus adjustment is made by moving the objective revolver. The stroke length is 7 mm.

Power switch

The illuminated power switch permits to switch the power supply of the microscope on and off. Due to the illuminated switch, the operating state of the microscope can easily be determined even in darkened rooms.

Fuse holder with voltage selector module

The fuse holder is equipped with two fuses and a voltage selector module. Depending on the country of use, the voltage needs to be set to 100 V, 115 V or 230 V (tolerance $\pm 10\%$).

Incident-light fluorescence facility*

The microscope version with incident-light fluorescence facility contains the integrated fluorescence axis and the lamp holder for mounting of a lamp housing.

Fluorescence filter block slide*

The fluorescence filter block slide holds up to 3 fluorescence filter blocks. The filter block slides can be moved to three alternative switching positions.

One position of the slide can also be used as a bright-field position by leaving it empty (no filter block inserted).

Lamps and lamp housings* for the incident-light fluorescence facility

For incident-light fluorescence, an additional illumination is required. On the DM IL, all lamp housings of the 106 and 107 ranges can be used. Depending on the version used, the controls for the lamp housings are located on the right or left side.

The LH 106 and LH 107 lamp housings are only used with a 12 V/100 W halogen lamp which can be centred in the X/Y plane. Both lamp housings are not equipped with a reflector but have diffusing screens, heat-absorbing filters and focusable aspheric collectors.

The LH 106z lamp housing is similar to the LH 106 lamp housing but has a reflector which can be centred and focused. It also contains a 4- or 6-lens collector (quartz collector on request).

The following lamps can be used with their special sockets:

Halogen lamp, 12 V/100 W

Hg ultra-high pressure lamp 50 W, AC

Hg ultra-high pressure lamp 100 W, DC, without igniter

Hg ultra-high pressure lamp 100 W, DC, with igniter ¹⁾

Xe ultra-high pressure lamp 75 W, DC, with igniter ¹⁾

¹⁾ This requires to raise the microscope stand with a base to increase the necessary clearance.

Power unit*

An external power unit is used to regulate the lamp for incident-light fluorescence and the corresponding lamp housings.

Modulation slide or phase contrast slide*

The modulation slide or phase contrast slide is part of a contrasting method, either Integrated Modulation Contrast (IMC) or phase contrast.

The same slides are used for the S 55 and S 90 condensers, however, with different phase rings.

If no phase slide or modulation slide is used, it is also possible to insert an empty slide into the corresponding holder on the condenser.


IMC modulator*

For the Leica Integrated Modulation Contrast technique, the IMC modulator is provided in the microscope stand. (→ p. 47 Operation of IMC). In the standard version, all DM IL microscopes are equipped with an empty slide.

Rear of the microscope

Fig. 4

- 1 6 V/35 W connection
- 2 Mains connection
- 3 Fuser holder with 2 mains fuses
- 4 Equipotential bonding point
- 5 Logo for equipotential bonding point
- 6 Lamp holder (for fluorescence version)
- 7 Safety note
- 8 Nameplate

Installation site

The microscope should be operated in a dust-free room which is also free from oil, chemical fumes, and excessive humidity. Major temperature fluctuations, direct sunlight, and vibration should be avoided because they can affect measurements and photomicrography.

Ambient conditions:

| | |
|-------------------|--------------------|
| Temperature | 10–36 °C |
| Relative humidity | 0–80 % up to 30 °C |

In warm and humid climates, microscopes require special care to prevent fungal growth. For additional instructions, please refer to chapters “Maintenance” and “Storage”.


Attention!

Lamp housings* and power units* must be located at least 10 cm away from the wall and combustible objects.

Unpacking

Please compare the delivery carefully with the packing note, delivery note or invoice. We strongly recommend that you keep a copy of these documents with the manual so that you have the proper information on the time and scope of delivery at a later time when ordering more equipment or when the microscope is serviced. Please make sure that no small parts are left in the packing material. Some of our packing material has symbols that indicate that it can be recycled in an environmental-friendly manner.

Note

Keep the packing material for the purpose of storage and transport of the microscope and its accessories.

Attention!

Avoid touching the surface of the optics lens. If finger prints are left on the glass surfaces, use a soft leather or linen cloth to remove them. Even minute traces of perspiration from fingers can corrode the surfaces of optical instruments in a very short time. For additional instructions, please refer to chapters “Maintenance” and “Cleaning”.

Attention!

On no account should you connect the power unit* and peripherals* to the mains yet!

The following parts can be included in the delivery:


- Leica DM IL microscope (transmitted-light or incident-light fluorescence)
- Illumination and condenser carrier
- Tube
- Eyepieces
- Objectives
- Condenser
- Protective cover
- Power cord
- Instruction manual

Optional components:

- IL/L tube adapter
- Phase slide
- Adjustment telescope
- IMC slide
- IMC module
- Filters for transmitted light
- Filter slide for fluorescence blocks
- Fluorescence blocks
- Lamp housing
- Replacement halogen lamp
- Ultra-high pressure mercury lamp
- External power unit
- c-mount video adapter
- Camera
- Speciment stage accessories
- Additional components from the DM L range, such as tubes, drawing facility, multi-discussion facility, magnification changer, ergo module

Fig. 5 DM IL microscope with drawing facility and ergo photo tube

1 Ergo photo tube from DML range, **2** IL/L tube adapter, **3** Drawing facility, **4** Specimen stage with accessories

Setting up

- As a first step, remove all components from the shipping and packing material.
- Place the DM IL basic microscope stand in a sufficiently clear space on a workbench.
- Make sure that all four feet have been mounted to the bottom of the microscope stand.
- If the delivery includes a stabiliser plate, mount this to the bottom of the microscope stand with two screws so that the two front feet fit into the recesses (Fig. 7). Tighten the screws and place the microscope back on the table in the upright position.

Attention!

On no account should you connect the microscope to the mains yet!

Attaching the condensers

- Screw the S 90 (8.1) or S 55 (8.2) condenser into the condenser holder (9.2) of the transmitted-light illumination carrier from below.

Fig. 6 Microscope with transmitted-light illumination column
1 Screw for condensor collision protection, 2 Feet


Fig. 7 Stabiliser plate



Fig. 8
1 S 90 condenser, 2 S 55 condenser

Positioning the transmitted-light illumination carrier

- Position the transmitted-light illumination carrier into the column from above while pressing on the notched lever for the condenser level adjustment (9.5).
- Position the transmitted-light illumination carrier (9.3) on the transmitted-light illumination column (9.4) depending on the condenser used (S 55 or S 90) and release the notched lever. The markings (9.6) are referred to a liquid level of 15 mm. For microscopes with dual markings, the lower line indicates a liquid level of 15 mm and the upper line indicates a liquid level of 50 mm.
- Make sure that the transmitted-light illumination carrier is locked in place.


Fig. 9 Transmitted-light illumination unit with condenser
1 Condenser (S 90), 2 Condenser holder, 3 Transmitted-light illumination carrier, 4 Transmitted-light illumination column, 5 Notched lever for condenser level adjustment, 6 Markings

Note

A screw (6.1) on the transmitted-light illumination column prevents that the condenser collides with the specimen stage.

Repositioning the transmitted-light illumination carrier with the specimen stage

Note

For the heating stage* only one position of the transmitted-light illumination carrier is possible. This position is determined by the existing holes in the stage.

Repositioning of the transmitted-light illumination carrier is only possible with the DM ILB or DM ILT tubes.

The specimen stage can be made accessible from three sides by repositioning it by 180°.

Attention!

Do not loosen the screws (10.1) below the stage on the transmitted-light illumination carrier because this would displace the optical axis.

- Unscrew the screws (11.1) with a 3 mm fixed spanner and remove the screws.
- Rotate the stage with the transmitted-light illumination unit by 180°.
- Position the stage with the transmitted-light illumination unit. The stage should lock in the guide pins (10.2) on the microscope.
- Insert the screws (12.1) and tighten them.
- Mount the connecting cable in the plastic guide below the stage (10.3).

Fig. 10

1 Clamping screws for optical axis, 2 Guide pins, 3 Plastic guide


Fig. 11

1 Screws for repositioning the transmitted-light illumination carrier



Fig. 12

1 Screws for repositioning the transmitted-light illumination carrier, 2 Connecting cable

Electrical connection of the transmitted-light illumination carrier

- Use the connecting cable (12.2) to connect the transmitted-light illumination with the integrated power supply via the power socket (13.1) on the rear of the instrument.

Positioning the tubes

The DM ILB tube (Fig. 14 Binocular tube) or the DMILT tube (Fig. 15 Binocular photo tube) is included in the standard delivery. For additional instructions, refer to chapter "Specifications; Performance data".

DM ILB and DM ILT tube

- Loosen the clamp (16.1) use a 3 mm fixed spanner.
- Insert the tube (16.3) into the tube holder (16.2).
- Retighten the clamp screw.
- To set a new viewing position, loosen the clamp screw (16.1) and retighten it after rotating the tube to the desired position.

Tube from the DM L range

In place of the standard DMIL tubes, you can also adapt one of the following tubes from the DM L range:

- HC LB 0/3/4 and HC LBP 0/3/4 binocular tube
- HC L1T 4/5/7 and HC L1TP 4/5/7 trinocular tube
- HC L3TP 4/5/7 trinocular tube with 3 switching positions
- HC LVB 0/4/4 ergo tube, binocular
- HC L1VT 0/4/4 ergo photo tube, trinocular

Fig. 13 Rear of microscope
1 Socket for connecting cable

Fig. 14 Binocular tube

Fig. 15 Binocular photo tube

Proceed as follows:

- Unscrew the lock screw on the tube changer (16.1) using a 3 mm fixed spanner.
- First place the IL/L tube adapter (17.1) into the tube holder of the microscope.
- Retighten the lock screw.
- Unscrew the lock screw on the IL/L tube adapter (17.2).
- Place a tube into the tube holder of the tube adapter.
- Retighten the lock screw.
- To set a new viewing position, loosen the lock screw (17.2), rotate the tube as desired and retighten the lock screw.

Note

For mounting a drawing facility*, a multi-discussion tube*, a magnification changer* or an ergo module*, please refer to chapter "Setting up the options".

Fig. 16

1 Lock screw, 2 Tube holder, 3 Tube


Fig. 17

1 IL/L tube adapter, 2 Lock screw

Fig. 18 HC L1T tube

Fig. 19 HC L1VT tube

Positioning the eyepieces

The eyepieces are positioned into the eyepiece connection pieces.

For the DM ILB and DM ILT tubes, use only the following eyepieces:

10x/18 Br (M),
10x/20 Br (M) or
15x/14 Br

When using a tube from the DM L range, use the corresponding eyepieces:

HC PLAN 10x/20 (M)
HC PLAN 12.5x/16 M
Widefield 16x/14 Br (M)⁺⁾
Widefield 25x/9.5 Br (M)⁺⁾

⁺⁾ also requires a spacer ring

Br) eyepiece for wearer of spectacles

For information on diameter, visible object surface and overall magnification of the microscope, refer to chapter "Specifications; Performance data".

Positioning the graticules*

Graticules can only be retrofitted by the user for the abovementioned HC PLAN eyepieces. For the abovementioned eyepieces for the DM ILB and DM ILT, the graticules need to be retrofitted at the factory.

Basically, graticules may only be used for eyepieces with adjustable eyelens = Type M.

Important:

Be extremely careful to avoid dust and fingermarks, since this will be visible in the field of view.

The graticule diameter is always 26 mm for HC PLAN eyepieces and 19 mm for the standard eyepieces for the DM ILB and DM ILT tubes.

Eyepieces 10x/18 M only:

- Unscrew the sleeve from the lower part of the eyepiece.
- Insert the graticule with the coated side pointing downwards (in the direction of the objective) so that any lettering appears the right way round when observed in the viewing direction later.
- Screw the sleeve back in.

Eyepieces HC PLAN 10x/20 M and HC PLAN 12.5x/16 M only:

- Unscrew the retainer sleeve from the lower part of the eyepiece.
- Insert the graticule with the coated side pointing downwards (in the direction of the objective) so that any lettering appears the right way round when observed in the viewing direction later.
- Screw the retainer sleeve back in.

Positioning the photo eyepieces*

The HC PLAN viewing eyepieces are designed for direct visual observation. Special eyepieces with an insertion diameter of **27** mm and the engraved identification **HC...PHOTO** are used for adapting micro-photographic equipment with fixed magnification, e.g. DMLD and MPS systems, as well as for special TV adaptation systems (use the proper adapter!).

Additional information on the adaptation of photo and video equipment → separate instructions.

Attaching and removing the objectives

- Remove the screw cap on the objective threads.
- Screw the objectives into the nosepiece opening that an incremental change of magnification is possible (e.g. in the sequence 4, 10, 20, 40).
- If objective threads are left unused, cover these with screw caps to protect the microscope optics against dust.

Positioning the filter


- Place the filter (Fig. 20) into the filter holder (21.1) on the transmitted-light illumination carrier.

Fig. 20 Filter

Fig. 22

- 1 Holes for attachable specimen controller, 2 Holes for specimen brackets, 3 Holes for stage mounting

Attaching the specimen controller

- Attach the specimen controller to the right or left side of the stage to accommodate brackets for different culture vessels (22.1).
- Mount the specimen controller with a 3 mm fixed spanner.

Fig. 21

- 1 Filter holder


Fig. 23 Specimen controller with holder mounting

Self-adhesive scales to determine the coordinate adjustment have been included with the brackets.

- Attach these in the recesses on the specimen controller.

Attaching the stage expansion

- Attach the stage expansions to the right or left side or to both sides of the stage (22.3).
- Mount the stage expansion with a 3 mm fixed spanner.

Inserting the specimen brackets

- Insert the specimen brackets into the four holes provided in the stage opening for this purpose (22.2).

Selecting the mains voltage and connecting the microscope to the mains supply

Attention!

Do not plug in the mains plug of the microscope and of the power unit* before all the options are assembled.


The mains plug may only be inserted into a grounded socket. The protection must not be jeopardised by using an extension cable without ground conductor.

Fig. 24 Specimen stage with expansions attached

Fig. 25

1 Fuse holder with voltage selector module, **2** Latch, **3** Voltage selector module, **4** Mains voltage, **5** Mains connection

- Check the voltage setting on the rear panel of the unit. Depending on the country of use, this setting can be 100 V, 115 V or 230 V. The voltage setting can be corrected as follows:
 - Use a screwdriver to press on the latch (25.2) and remove the fuse holder (25.1).
 - Pull out the voltage selector module (25.3).
 - Insert the voltage selector module into the holder so that the number which corresponds to the desired mains voltage (25.4) appears on the outside (upside down).
 - Insert the fuse holder (2.1) until the latch clicks audibly in place.
- If you purchased additional options with the microscope, install these options first (see next chapter).

Attention!

For external lamp power units, always perform the mains voltage setting as described in the manual which is supplied separately, or use a power transformer.

- Then connect the microscope with the power cord (2.4) and connect it to the mains supply.

Attention!

Observe the safety information on pages 7 and 8!

Setting up the options

Note

These assembly steps are not required if no additional accessories have been purchased with the microscope.

Assembling the fluorescence filter blocks*

Note

For microscopes with integrated incident-light fluorescence facility only.

The filter block slide (Fig. 26) accommodates up to three fluorescence filter blocks.

- For assembly of the filter blocks, remove the slide cover (26.5).
- Insert the filter blocks (26.1) into the dovetail holder (26.2), with the engraved lettering pointing upward. Make sure that the filter blocks click in place.
- Reinstall the slide cover.
- Check that the slide cover (26.5) is properly installed.
- Mark the filter positions using the enclosed adhesive labels.

Positioning the filter block slide*

- Hold the filter block slide so that the warning label is positioned is located at the front left side and insert it into the dovetail holder (26.2) on the left side of the microscope stand.
- The filter block slide can now be moved between the three switching positions.

Attention!

The filter block slide is not protected against inadvertent sliding out.

Attention!

When using an incident-light dichroic mirror or polarisation dichroic mirror in conjunction with fluorescence filter blocks, there is a risk of glare in case of inadvertent switching between positions!

Positioning the phase contrast light rings

The phase contrast light rings are inserted into the phase contrast slide (Fig. 27). The slide has two holders which can be centred and a centre position which is the bright-field position (BF).

- Place the light rings into the centring holders (27.2) of the slide with the disc pointing downwards.
- Press on the edge of the light rings until they click in place. Avoid any pressure on the centre of the disc because this can break the webs.

Different sets of light rings are available for the S 55 or S 90 condensors.

Positioning the phase contrast slide on the transmitted-light illumination carrier*

! Attention!

The condenser slide is not protected against inadvertent sliding out.

- Remove the blank slide if present.
- Hold the condenser slide so that the lettering points forward. The catches (27.4) are located on the **upper** long side of the slide and point toward the centre of the specimen stage.
- Insert the condenser slide into the transmitted-light illumination carrier from the side (Fig. 28). The keyways should click in place when the slide is inserted.

Fig. 26

1 Filter block, 2 Dovetail holder, 3 Slide cover, 4 Lower part of filter block slide


Fig. 27

1 Slide for light rings, 2 Centring holder for light ring, 3 Bright-field position, 4 Catch



Fig. 28 Transmitted-light illumination carrier with holder for condenser slide

Positioning the IMC slit diaphragm slide* on the transmitted-light illumination carrier

! **Attention!**

The condenser slide is not protected against inadvertent sliding out.

- Remove the blank slide if present.
- Hold the condenser slide so that the lettering "Top left" is located on the left upper side (29.1) and the other lettering points forward. The catches (29.2) are located on the **upper** long side of the slide and point toward the centre of the specimen stage.
- Insert the condenser slide into the transmitted-light illumination carrier from the side (30.1). The keyways should click in place when the slide is inserted.

Positioning the IMC modulator*

- Remove the blank slide if present.
- Insert the IMC modulator so that the lettering (29.3) points forward.
- Lock the slide in position BF (lettering BF is visible) or in position IMC (lettering IMC is visible).

Fig. 29 IMC slit diaphragm slide and IMC modulator
1 Lettering "Top left" on condenser slide, 2 Catches, 3 Lettering on IMC modulator



Fig. 30
1 Holder for IMC modulator

Assembling the 106* or 107* lamp housings

Note

For microscope with integrated incident-light fluorescence facility only.

Attention!

Prior to making any assembly work, always disconnect the power supply at the external transformer outlet and the microscope outlet!

Replacing or changing the halogen lamp

This procedure is only required if the lamp is not pre-installed.

- Unscrew the screw on the lid (31.1, 32.1) with a crosstip screwdriver and remove the lid (Fig. 33).
- Move the collector to the front.

Note

This step is not required for the 107/2 lamp housing.

Attention!

Leave the protective covering on the lamp until the lamp is in its holder! Avoid making finger prints or wipe them off immediately.

Fig. 31 Series LH 106 lamp housing

1 Screw for opening of lamp housing, 2, 3 X and Y centring of lamp*, 4 Focusing of collector, 5, 7 Mounting screws, 6 Filter holder (intermediate component) for filter Ø 50 mm



Fig. 32/33 Opening the LH 107 lamp housing

1 Screw for opening of lamp housing

- Place a new 12 V/100 W halogen lamp straight into the lamp holder (34.1).
- Move the collector back to its original position.
- Replace the lid and fix it with the screw (31.1, 32.1).
- Connect the lamp housing to the power unit:

Assembling the 106 z* lamp housing with halogen lamp

Note

For microscope with integrated incident-light fluorescence facility only.

Attention!

Prior to making any assembly work, always disconnect the power supply at the external transformer outlet and the microscope outlet!

- Unscrew the lock screws (36.4 and 36.9) with a crosstip screwdriver.
- Slightly pull out the cut-out plug (36.11) from the socket and flip up the lid (36.1).

Attention!

Leave the protective covering on the lamp until the lamp is in its holder! Avoid making finger prints or wipe them off immediately.

Fig. 34 LH 106 lamp housing, opened

1 Lamp holder with 12 V/100 W halogen lamp, 2 Collector, 3 Diffuser



Fig. 35 107/2 lamp housing, opened

1 Collector, 2 Lamp holder with 12 V/100 W halogen lamp

- Unscrew the lock screws (36.10) on the lamp holder and remove the lamp holder (Fig. 37).
- Insert a new 12 V/100 W halogen lamp into the lamp holder.
- Insert the lamp holder and fix it with the screws (36.10).
- Plug the cut-out plug into the socket (36.11).
- Flip the lid back down and tighten the screws (36.4 and 36.9) on the lid.
- Mount the lamp housing and fix it on the microscope with the lock screw.
- Connect the lamp housing to the power unit.


Connecting to the 12 V/100 W power unit*

The Leica 12 V/100 W power unit is an autoranging power supply. It supplies halogen lamps up to 12 V.

Be careful when unpacking the power unit. Do not operate the unit in the presence of strong vibration, large temperature fluctuations, high humidity and direct sunlight. Position the unit so that the ventilation openings are clear. Keep a safety distance of 10 cm from walls or other equipment.

Fig. 36 106 z lamp housing, opened

1 Lid, flipped up, **2** Collector, **3** 12 V/100 W halogen lamp or gas discharge lamp (see Fig. 40), **4, 9** Lid mounting, **5** Reflector, **6, 8** Adjusting screws for x/y centring of reflector, **7** Focusing of reflector, **10** Mounting screws for lamp holder, **11** Socket for cut-out plug

A rotary knob is located on the front panel of the power unit which is used to set the desired voltage. The setting range is between approx. 2.5 V and approx. 12 V. The 10.5 V setting has a marking and corresponds to a colour temperature of 3200 K. This position has been preset at the factory.

The On/Off switch and the ports for the mains plug and the lamp are located on the rear panel of the unit. The "Control" port is currently not used.

Connect the lamp housing to the corresponding port. When using a lamp housing with shielded cable (e.g. the 107/2 lamp housing), screw the shield connection to the equipotential bonding point. Connect the power unit to the mains. It is not required to preselect the mains voltage. Switch on the unit with the On/Off switch (Position 1).

Assembling the 106 z* lamp housing with Hg and Xe lamps

Note

For microscope with integrated incident-light fluorescence facility only.

The following gas discharge lamps (with different lamp holders) can be used in addition to halogen lamps (Fig. 40):

Hg ultra-high pressure lamp 50 W, AC

Xe ultra-high pressure lamp 75 W, DC, stabilised

Hg ultra-high pressure lamp 100 W, DC, stabilised

Hg ultra-high pressure lamp 100 W, DC, stabilised, with igniter

Fig. 37 Front panel of 12 V/100 W Leica

Fig. 38 Rear panel of 12 V/100 W Leica
1 Equipotential bonding point

Attention!

- Prior to making any assembly work, always disconnect the power supply at the external transformer outlet and the microscope outlet!
 - Allow the lamp housing to cool prior to opening it (at least 15 min.), explosion hazard!
 - Never touch the glass parts of the burner with your hands. If required, carefully remove finger prints and dust (use alcohol if necessary).
 - Adjust lamps immediately after ignition.
 - Avoid frequent switching on and off because this could affect the service life (→ p. 76) and stability. Hot Hg lamps will ignite only after cooling off. It is recommended to run in new burners for a couple of hours without interruption.
 - Make sure that the lamp housing is sufficiently ventilated. Never block air slots with paper etc., fire hazard!
 - It is good practice to record the hours of use and to compare them to the manufacturer's specifications.
Replace discoloured, worn burners in due time.
 - We must refuse any liability for damage resulting from a possible explosion of the lamp.
-
- If required, disconnect the mains plug of the power unit and the microscope.
 - Open the 106z lamp housing by slightly loosening the screws (36.4), pulling out the cut-out plug from the socket (36.11) and flipping up the lid of the lamp housing.
 - Unscrew the safety screws (36.10) and pull out the lamp holder (Fig. 39, 40).

- Insert the burner as follows while strictly observing the above safety instructions:
 - If a plastic protective sleeve is present, leave it in place for the time being.
 - Insert the burner so that the lettering is **upright** after installation. For Hg 50, Hg 100, and Xe 75, the different height of the metal base ensures installation at the proper height.
 - Align any existing glass fused nipple (40.2) by rotating the burner so that the nipple is orientated **to the side** and away from the beam path.
 - Insert the upper pin of the burner between the clamps of the flexible power supply and fix it with the screw (40.1) .
 - Slightly unscrew the stud (40.4) in the holder.
 - Insert the burner into the lower end of the metal base and retighten the stud.
 - Remove the protective sleeve of the burner if it is still present.

Fig. 39 12 V/100 W lamp holder

- Place the lamp holder with the burner inserted into the lamp housing and tighten the screws (36.10).
- Close the lid of the lamp housing.
When closing the lamp housing, make sure that the pins of the cut-out plug engage in the sockets provided for this purpose.
- Retighten the screws of the lid.
- Push the cut-out plug fully in.
- Attach the lamp housing and fasten it to the microscope with the lock screw.
- Connect the lamp housing to the power unit (compare mains voltage!):

The Hg 50 W lamp is properly installed if:

1. The type is stamped on the lower socket of the lamp. The stamped lettering should be visible, i. e. **not** upside down.
2. The upper base is marked "UP".

Note:

An incorrect installation will reduce the lamp brightness to about 60 % and will considerably limit the useful life of the lamp.

Attention:

Make sure that the markings on the lamp base and on the power unit are the same. For example, if the lamp base is marked L1, L1 must also be set on the power unit to make full use of the lamp and not to shorten its life.


Important:

Make sure to dispose of worn burners in an environmental-friendly manner.

Fig. 40 Lamp holders for gas discharge lamps

1 Upper clamp, **2** Fused nipple of burner, **3** Lower clamp, **4, 6** Mounting holes for holder, **5** Sockets for cut-out plug, **7** Protective sleeve

Adaptation of imaging systems on binocular photo tubes*

The binocular photo tubes provided for the Leica DM IL can be used to adapt an imaging system, e. g. video camera, reflex camera or automatic photo system (e. g. Leica MPS 48/52).

Microphotography

Requirements for the adaptation of microphotographic equipment include a trinocular tube, a HC PHOTO eyepiece adapter tube and HC PHOTO eyepieces with 27 mm insertion diameter. If the microphotographic equipment is not provided with a special viewing port to limit the format, it is necessary to use HC PLAN M eyepieces, i. e. with focusable eyelens and photo graticule inserted, in the binocular viewing port. For additional details, please refer to the manual supplied with the photo equipment.

TV adaptation

Different adapters are available for connecting TV cameras with c-mount and B-mount objective thread. The adapters listed in the table below can be used on all trinocular photo tubes, although some tubes require an additional photo adapter tube. The picture area on the monitor depends on the adapter used and on the chip size of the camera.

Calculation of the magnification on the monitor

The magnification V_{TV} on the monitor can be calculated using the equation below or measured with an object micrometer and a cm scale.

$$V_{TV} = \frac{\text{Objective magnification} \times \text{Magnification changer* factor} \times \text{TV adapter magnification} \times \text{Monitor diameter}}{\text{Chip diameter of camera}}$$

Recorded picture diagonal in mm for
1-inch camera 2/3-inch camera 1/2-inch camera 1/3-inch camera

Without zoom magnification:

| | | | | |
|---|----|------|------|------|
| c-mount adapter 1x HC | 16 | 11 | 8 | 6 |
| c-mount adapter 0.63x HC ⁺⁾ | – | 17.5 | 12.7 | 9.5 |
| c-mount adapter 0.5x HC | – | – | 16 | 12 |
| c-mount adapter 0.35x HC | – | – | – | 17.1 |
| c-mount adapter 4x HC ⁺⁾ | 4 | 2.8 | 2 | 1.5 |

With zoom magnification (Vario TV adapter):

| | | | | |
|---------------------------------------|---|---|-----------------------|----------|
| c-mount, 0.32 – 1.6x HC | – | – | 19 ⁺⁺⁾ – 5 | 18 – 3.8 |
| B-mount, 0.5 – 2.4x HC | – | – | 16 – 3.3 | – |
| B-mount, 0.5 – 2.4x HC ⁺⁾ | – | – | – | 12 – 2.5 |

⁺⁾ in preparation


⁺⁺⁾ zoom factor 0.42 x and higher only!

The procedure for positioning the drawing facility, the multi-discussion facility, the magnification changer, and the ergo module is basically the same. The drawing facility or the ergo module can either be mounted directly on the basic microscope stand (in conjunction with the DM ILB or DM ILT tubes) or on the DM IL/L tube adapter (in conjunction with L-tubes).

Positioning the drawing facility*

- Unscrew the lock screw (41.1) on the microscope stand using a 3 mm fixed spanner.
- When using the DM IL/L tube adapter:
 - Attach the IL/L tube adapter (41.2).
 - Retighten the lock screw (41.1).
 - Unscrew the lock screw (41.4) on the tube adapter.
- Insert the drawing facility (41.3) into the tube holder of the basic microscope or of the IL/L tube adapter.
- Retighten the lock screw (41.4).
- Unscrew the lock screw on the drawing facility.
- Position the tube.
- Retighten the lock screw on the drawing facility.

Fig. 41 Microscope with drawing facility
 1 Lock screw, 2 IL/L tube adapter, 3 Drawing facility, 4 Lock screw

Positioning the multi-discussion facility*


Note

When using the multi-discussion facility, the stabiliser plate should also be mounted (→ p. 17).

- Unscrew the lock screw (41.1) on the microscope stand using a 3 mm fixed spanner.
- Attach the IL/L tube adapter (41.2).
- Retighten the lock screw (41.1).
- Unscrew the lock screw (41.4) on the tube adapter.
- Insert the multi-discussion facility (42.1) into the tube holder of the adapter.
- Retighten the lock screw (41.4).
- Unscrew the lock screw on the multi-discussion facility.
- Position the tube.
- Retighten the lock screw on the multi-discussion facility.

Fig. 42 Microscope with multi-discussion facility
 1 Multi-discussion facility

Positioning the magnification changer* (not shown)

- Unscrew the lock screw (41.1) on the microscope stand with 3 mm fixed spanner.
- Attach the tube adapter IL/L (41.2).
- Retighten the lock screw (41.1).
- Unscrew the lock screw (41.4) on the tube adapter.
- Insert the magnification changer into the tube holder of the adapter.
- Retighten the lock screw (41.4).
- Unscrew the lock screw on the magnification changer.
- Position the tube.
- Retighten the lock screw on the magnification changer.

Positioning the ergo module* (not shown)

- Unscrew the lock screw (41.1) on the microscope stand with 3 mm fixed spanner.
- When using the DM IL/L tube adapter:
 - Attach the IL/L tube adapter (41.2).
 - Retighten the lock screw (41.1).
 - Unscrew the lock screw (41.4) on the tube adapter.
- Insert the ergo module into the tube holder of the basic microscope or of the IL/L tube adapter.
- Retighten the lock screw (41.4).
- Unscrew the lock screw on the ergo module.
- Position the tube.
- Retighten the lock screw on the ergo module.

Operation

Attention!

Be especially careful when making examinations which involve the use of acids or other aggressive chemicals. Avoid direct contact of these substances with optical and mechanical parts.

Basic setting for transmitted light

Switching on the 6 V/35 W halogen lamp

- Switch on the 6 V/35 W halogen lamp with the power switch (43.2) (position I = On).
- Adjust the brightness with the rotary knob (43.1).

Adjustment specimen

For the initial adjustment of the microscope it is recommended to use a specimen that contains areas of high and low contrast.

For incident-light fluorescence of transparent objects, adjust in transmitted light first.

Fig. 43

1 Brightness adjustment knob, **2** Power switch, **3** Focusing mechanism

Focusing on the object

- Set the desired objective. For this purpose, lower the objective nosepiece first. Use the black knurled grip on the nosepiece to move the objective into the light path. Make sure that the nosepiece clicks audibly in place.
- Use the coarse and fine adjustment to focus on the object; this will vary the height of the objective nosepiece while the stage level will remain the same. The overall travel is 7 mm. The focusing range (in air) is from 1.0 mm below to 6 mm above the stage surface.

Attention!

Depending on the objective used, the objective nosepiece **must** be lowered before changing the objective position. Otherwise the objective may collide with the stage.

Adjusting tubes and eyepieces

If you wear glasses, you should remove or fold back the glare protection on the eyepieces but make sure to put it on if you are not wearing glasses.

- Set your interpupillary distance by pulling the eyepiece tubes apart or pulling them closer together until you see a single superimposed image instead of a double image.
- Note your personal interpupillary distance.
- Additional procedure for ergo tubes: Set the viewing angle ($0^\circ - 35^\circ$) by tilting the binocular viewing port. To avoid fatigue symptoms, vary the viewing angle from time to time.
- Close any unused tube exits as stray light may otherwise disturb viewing.

DM ILB binocular tube

Only for eyepieces with graticule* inserted:

- Largely defocus the object or remove it from the beam path.
- Focus the graticule with the eye relaxed by adjusting the eyelens. (The eye relaxes most if you look temporarily at a distant object outside the frame.)
- Focus the object only through the eyepiece with the graticule.
- Then close the eye and focus the object by adjusting the second eyepiece only.

Only when no graticule is inserted in both eyepieces:

- Largely defocus the object or remove it from the beam path.
- Adjust the eyelens so as to focus on the boundary of the field of view. When adjusting the eyelens you will see a bright line around the body of the eyepiece. This line shows the correct position of the eye-lens for people with normal eyesight and for wearers of glasses who look through the microscope with corrective glasses.

Remove glasses with bifocal or progressive lenses before looking through the microscope.

- Focus the object through the eyepieces.

Only if one eyepiece has no adjustable eyelens:

- Focus the object through this eyepiece first (close your other eye).
- Focus the image again by adjusting the eyelens of the second eyepiece.

Correction for defective vision

- Look through the right eyepiece with your right eye and use the fine adjustment to focus on the specimen.
- Then look onto the same spot of the specimen with your left eye and rotate the left eyepiece connection piece until the point on the object is sharply focused. Do not use the fine adjustment to achieve this.
- If eyepieces with adjustable eyelenses are used, do not correct for defective vision by adjusting an eyepiece tube but by adjusting the eyelens of the eyepiece.

DM ILT trinocular tube

- Switch the beam splitter to visual observation by moving the rod. The meaning of the switching positions is marked with symbols on the lateral face of the tube.

Rod retracted = visual

Rod inserted = photo

- Adjustment of the eyepieces is made in the same manner as for the binocular tube.
- Correct for defective vision by adjusting the eyelens of the eyepiece.

Fig. 44 DM ILB tube

Fig. 45 DM ILT tube

Operation of objectives

Immersion objectives

OIL: Use optical immersion oil in compliance with DIN/ISO only.

Cleaning → p. 57, Labelling → p. 67 pp.

Attention:

Observe the safety instructions for immersion oil!

W: Water immersion. The special water immersion objectives with a ceramic front can be used for all aqueous solutions.

IMM: Universal objective for water, glycerine, oil.

Colour codes of objectives → p. 69

Locking of objectives

Some immersion objectives (identified by a knurled grip) can be virtually shortened (locked). This stops any remaining drops of immersion liquid from inadvertently wetting objectives and other objects.

- Push in the front part by about 2 mm in the direction of the nosepiece.
- Lock the objective by rotating it slightly.

Attention:

The locking mechanism must always be released before the immersion objective is used again, since otherwise the spring mechanism protecting the specimen and the objective is ineffective and the other objectives are no longer parfocal with the immersion objective.

CORR objectives

These are special objectives which can be adjusted to different coverglass thicknesses.

- Adjust the correction mount coarsely to an average or estimated value by using the knurled grip.
- Focus the preparation.
- Readjust the correction mount until you obtain an optimum contrast, which may require refocusing with the fine adjustment.

Operation of transmitted light

Bright-field illumination

Illumination methods which show object-free areas of the specimen as the brightest parts of the image are called bright-field. Bright-field illumination requires absorbing object structures, i. e. staining of the specimen is useful in most cases. Alternatives are optical contrast methods, such as phase contrast or modulation contrast.

Adjusting the condenser

To aid in proper level adjustment of the S 90 und S 55 condensers, markings (46.1) are located on the column. These markings indicate a liquid level of 15 mm. For microscope stands with dual markings, the lower line indicates a liquid level of 15 mm while the upper line indicates a liquid level of 50 mm.

- Press the catch lever (46.2) and adjust the incident-light illumination carrier until the upper edge of the carrier and the corresponding condenser level marking match.

Adjusting the aperture diaphragm

The aperture diaphragm (46.3) determines resolution, depth of field, and contrast of the microscope image. The best resolution is obtained when the apertures of the objective and of the condenser are roughly the same. Reducing the aperture diaphragm to be smaller than the objective aperture will reduce the resolution but enhance the contrast. A noticeable reduction in resolution occurs when the aperture diaphragm is reduced to less than 0.6x of the objective aperture and should be avoided where possible.

- Adjust the aperture diaphragm according to your subjective impression of the image.
- You may basically achieve a calibration by comparison with the apertures of different objectives.
- Visual comparison between the apertures of the objective and the condenser can be made as follows:
 - Remove the eyepiece from the eyepiece tube or use an focusing telescope and focus.
 - Close or open the aperture diaphragm until its image is just visible in the objective pupil (= brighter circle). This is considered the standard setting, i. e. condenser aperture = objective aperture.
 - Reattach the eyepiece.


For low-contrast objects, the aperture diaphragm can be closed further to highlight even faint specimen details. In polarised light microscopy, narrowing the aperture diaphragm usually results in brighter colours.

Attention:

The aperture diaphragm in the **illumination beam path** is **not** for setting the image brightness. Only the rotary brightness adjustment knob or the neutral density filter should be used for this purpose.

Fig. 46 Transmitted-light illumination unit with condenser
1 Markings, 2 Catch lever for condenser adjustment, 3 Aperture diaphragm

An aperture diaphragm in the **objective** is usually fully opened. The reduction in image brightness caused by stopping down results in:

- Better depth of field
- Lower coverglass sensitivity
- Suitability for darkfield
- Change in contrast

Possible errors

Wrong coverglass thickness or wrong objective.
Specimen has been placed on the stage with the coverglass upwards instead of downwards.

Aperture diaphragm opened too wide or closed.

Condenser at incorrect level.

Light ring inadvertently used.

IMC component inadvertently used.

Dirty optics.

Operation of phase contrast

Phase contrast is used to form high-contrast images of unstained specimens.

- Adjust the condenser level.
- Place the slide with the required light rings into the holder (47.1).
- Rotate the objective nosepiece to move the phase contrast objective (engraving PH) with the lowest magnification into the light path.
- Open the aperture diaphragm (47.4) marked "PH".
- Focus the specimen using coarse and fine adjustment. If it proves difficult to find the object plane: Temporarily narrow the aperture diaphragm or use a stained specimen. Set the condenser disc to the BF position or pull out the light ring slide. Reopen the aperture diaphragm.
- Use the light ring (e. g. **1**) that corresponds to the engraving on the objective (e. g. PH **1**).
- Centre the light ring as follows:
 - Remove one eyepiece from the tube.
 - Insert the focusing telescope.
 - Loosen the clamp ring on the focusing telescope and shift it until the light ring (bright) and the phase ring (dark) are in sharp focus.
 - If the light and the phase rings are greatly different in dimension, they need to be matched by varying the condenser level.
 - If the light ring is offset laterally against the phase ring, centre the light ring. To do so, rotate the centring key in the centring screws (47.5) until the phase ring covers the light ring.

Possible errors

Specimen: too thick, too thin, too brightly stained; identical refractive index of mounting medium and specimen so that no phase jump occurs.

Attention:

Wedge-shaped coverglass position which renders the centration of light and phase ring ineffective.

Wrong light ring, or light ring has been installed at wrong level.

Aperture diaphragm not opened.

Light ring not centred.

Wrong light ring slide.

IMC modulator in IMC position.

Condenser S 55 and condensor S 90 exchanged.

Fig. 47

1 Holder for light ring and light segment slide, **2** Transmitted-light illumination column, **3** Catch lever for condenser level adjustment, **4** Aperture diaphragm, **5** Centring screws for light rings, **6** Filter holder \varnothing 32 mm



Fig. 48

1 Light ring offset against phase ring: no phase contrast, **2** Phase ring fully covers the light ring: phase contrast

Operation of Integrated Modulation Contrast (IMC)

Integrated Modulation Contrast is a special form of oblique illumination based on the principle of Hoffman's modulation contrast.

With this method, a modulator is used to convert the phase gradients of an unstained object into amplitude differences.

The impression of a three-dimensional image is created, similar to a microscopic image with interference contrast. However, other than with interference contrast, the object can also be viewed through double-refracting plastic materials, such as Petri dishes.

Additional benefits of this imaging method include:

- High contrast
- High resolution
- Halo-free, variable-contrast relief image
- Long working distance of the condenser
- Simple assembly and adjustment
- Applicable for stained and unstained specimens

Important!

IMC is only possible in conjunction with the S 55 condenser.

Standard bright-field and phase contrast objectives can be used for IMC, which permits to cover the magnification range from 5x to 100x.

The following objectives are especially suited:

C PLAN 10x/0.22 AP 32.2

C PLAN L 20x/0.30 D

C PLAN L 40x/0.50 D

plus the corresponding phase contrast objectives.

All other objectives with pupil position D can also be used.

The following objectives with pupil position C may also be used with some restrictions:

N PLAN L 20x/0.40 Corr

N PLAN L 40x/0.55 Corr

PL FLUOTAR® L 63x/0.70 Corr

(Refer also to Possible errors).

The IMC requires the use of the IMC modulator (49.1) and the IMC slit-diaphragm slide (49.2).

Positioning the IMC modulator

- Remove the empty slide in the microscope stand if present.
- Position the IMC modulator so that the lettering points forward.
- Lock the slide in position IMC (lettering IMC visible).

The IMC modulator will be flush on either side.

Positioning the IMC slit-diaphragm slide on the transmitted-light illumination carrier

- Remove the empty slide or the phase contrast slide in the transmitted-light illumination carrier, if present.
- Hold the diaphragm holder so that the lettering "Top left" is located on the top left side and the other lettering points forward. The catches are located on the **upper** long side of the slide and point toward the centre of the specimen stage.
- Insert the diaphragm slide from the right side into the transmitted-light illumination carrier.
- Lock the slide in position IMC (lettering IMC visible).

Adjusting the slit-diaphragm

- Fully open the aperture diaphragm.
- Select a medium brightness since the line of light would otherwise appear too bright.
- Switch off any filters which may be switched on.
- Rotate the objective with the lowest magnification into the beam path, which will usually be the 10x objective.
- To adjust the slit width, move the slide of the IMC slit-diaphragm slide to the proper position for the objective, e. g. to the position marked 10x for the 10x objective.
- Remove one eyepiece and insert the focusing telescope.
- The line of light will appear as a bright line on the grey image of the modulator. Focus the line of light with the focusing telescope.
- Adjust the position of the line of light with the adjusting screws to the right of the IMC slit-diaphragm slide. A fitting Allen key is supplied for this purpose.


Attention:

Do not loosen the screws (49.4) on the slide!

- The line of light must be fully located on the grey field. With the 10x objective, the image of the modulator and of the line of light have almost the same size. Adjust the slit-diaphragm so that the boundary of the bright line of light is located near the darker edge.
- Move the other objectives with ascending magnification into the light path one after the other and check the position of the line of light. In the case of smaller deviations, find an intermediate position.
Always make sure that the objective magnification and the slide position on the IMC slit-diaphragm slide match.

Fig. 49 IMC components

1 IMC modulator, **2** IMC slit-diaphragm slide, **3** Catches, **4** Screws

Optimising the IMC (Fig. 49a):

When using the objective with the largest magnification, it may not be possible to adjust the light slit stop perfectly (in other words, the light slit cannot be positioned completely on the grey field, with the result that an offset is noticeable in either the white or the dark area). In this case, it is also possible to perform fine adjustments using the adjusting screw on the IMC modulator slide.

Use the Allen key to turn this screw to minimise the offset (to superimpose the grey area and the illumination slit). After that, swivel the 10x objective back in (followed by the other objectives) and adjust at the modulator again as described above. Repeating this operation several times should mean that offset no longer occurs. This adjustment generally need to be performed only once. Once the IMC is perfectly adjusted, remove the focusing telescope and replace the eyepiece.


Fig. 49a IMC-Components (new)

1 IMC-Modulator, **2** IMC slit-diaphragm slide, **3** Allen key, **4** Adjusting screw on the IMC slit-diaphragm slide, **5** Adjusting screw on the IMC-Modulator-slide

Possible errors

Poor image quality due to use of objectives without pupil position D.

Try improving the image quality as follows:

Reverse the IMC modulator (lettering to the rear). Adjust the slit so that sufficient coverage is achieved to avoid glare.

Short of optimum position of the slit-diaphragm.

The IMC modulator or the IMC slit-diaphragm slide are not locked in the IMC position.

Incorrect condenser level or wrong condenser (only S 55 condenser possible!).

Fluorescence filters are not disabled.

Operation of incident-light fluorescence

Note

For microscopes with integrated incident-light fluorescence facility only.

For the viewing of transparent objects using incident-light fluorescence, it is recommended to make an adjustment with transmitted-light first.


- Open the light stop by moving the lever (50.1).
 - = Light stop moved out of beam path
 - = Light stop moved into beam path
- Slide the filter block into the beam path (50.3).
- Position the specimen and focus. The field diaphragm is installed and pre-centred so that no adjustment is necessary.

If a reddish background of the specimen becomes visible with UV excitation, you can eliminate this effect by moving the BG 38 red attenuating filter (50.2) into the beam path.

- = Filter moved out of beam path
- = Filter moved into beam path

Fig. 50

1 Light stop, 2 BG 38 filter, 3 Filter block slide

Switching on and adjusting the 12 V/100 W halogen lamp in the 106* lamp housing

- Switch on the 12 V/100 W halogen lamp at the power unit.
- Open the light stop.
- Move the filter block into the beam path.
- Remove the objective in the beam path.
- Place a white sheet of paper on the specimen stage.
- Rotate the collector adjustment (51.4) until the lamp filament (Fig. 52) is clearly projected.
- Use a 3 mm fixed spanner to adjust the centring screws for level adjustment (51.2) and for horizontal adjustment (51.3) of the lamp until the lamp filament is in the centre of the light spot.

- Remove the paper.
- Position the specimen.
- Check with low objective magnification that the image is illuminated in a homogeneous manner.
- Adjust the collector (51.4) if necessary.

Switching on and adjusting the halogen, Xe and Hg lamps in the 106 z lamp housing*

- Switch on the lamp at the power unit.
- Open the light stop.
- Move the filter block into the beam path.
- Place a white sheet of paper on the specimen stage.
- Coarsely focus on the surface using a dry objective with low to medium magnification.

Fig. 51 106 lamp housing (with 12 V/100 W halogen lamp)

1 Screw for opening the lamp housing, 2, 3 X/Y centring of the lamp (compartment for storing a 3 mm allen key or screw-driver), 4 Collector focusing, 5, 7 Lock screw for mounting, 6 Filter holder (intermediate piece) for 50 mm filter
Items 3–4 do not apply to the 107 lamp housing



Fig. 52 106 lamp housing

Reflection of the lamp filament, greatly schematised: In reality, the reflection is extremely low in contrast, the bright overlap area is wider and less defined. With the 106 z lamp housing, the reflection is rotated by 90°.

- Use a pen to draw a mark in the centre of the bright surface.
- Remove the objective in the beam path.
- Rotate the collector adjustment (53.6) until the lamp filament or the discharge arc is clearly projected.
- Move the **reflection** of the lamp filament or the discharge arc to the side (54 a) by rotating the adjustment screws on the rear of the lamp housing (53.2 and 53.4).
- Focus the **direct image** of the lamp filament or the discharge arc and adjust it as follows:

For halogen lamp:


- Move the direct image to a position just below or above your centre mark (54b) or, especially for higher objective magnifications such as with the Xe lamp (54c), to the centre, i.e. superimposed.
- Move the reflection into the brighter circular area.
- Align the reflection symmetrically with the direct image (54c). Alternatively, you may also superimpose the two images as for the Hg and Xe lamps.

For mercury (Hg) and xenon lamps (Xe):

- Use the horizontal and vertical adjustment controls of the holder (53.5 and 53.1) to move the direct image into the centre of the brighter circular area.
- Move the reflection into the brighter circular area.
- Focus the reflection.
- Adjust the mirror until the reflection is superimposed to the direct image (54c).
- Remove the paper.
- Position the specimen.
- Check with low objective magnification that the image is illuminated in a homogeneous manner.
- Adjust the collector (53.6) if necessary.

Fig. 53 106 z lamp housing


1 Level adjustment of lamp, **2, 4** Level and lateral adjustment of reflection, **3** Mirror focusing, **5** Lateral adjustment of lamp, **6** Collector (focusing of lamp image), **7** Mounting screw

Attention!

Be careful not to project the reflection onto the electrodes as there is a risk of overheating leading to explosion. The two electrodes can just be seen in the extension of the symmetry plane of the discharge arc.

Fig. 54 Schematic view for 106 z lamp housing (in reality, lamp images are less defined)
a direct lamp image focused but out of centre
b direct lamp image in desired position
c indirect and direct lamp images in desired position

Possible errors

Weak fluorescence, weak image intensity:

Specimens improperly stored, too old or faded.

Rapid fading of specimens (e. g. for FITC).

Unspecific filter combination.

Numerical aperture of objectives too low.

Eyepiece magnification too high.

Spent lamp.

Room too bright.

Trinocular tube: incorrect beam splitter setting.

Secondary light due to reflection at condenser.

Low-contrast image due to:

Excitation bandwidth too wide.

Unspecific staining.

Fluorescing mounting medium.

Auto-fluorescence of objective or immersion oil.

Dirty glass surfaces.

Care and maintenance

Attention!

Disconnect the mains plug before cleaning and servicing!
Protect electrical components from humidity!

In warm and humid climates, microscopes require special care to prevent fungal growth. Clean the microscope after each use and make sure to keep the microscope immaculately clean.

Dust protection

Note

Protect the microscope and its accessories from dust by putting on the dust cover after each work session.

Cleaning

Attention!

Fibre and dust residues can cause disturbing background fluorescence during fluorescence microscopy.

Cleaning painted parts

Dust and loose particles of dirt can be removed with a soft brush or lint-free cotton cloth.

Obstinate dirt can be removed with commercially available aqueous solutions, benzine or alcohol.

To clean painted parts, use a linen or leather cloth moistened with any of these agents.

Attention!

Do not use acetone, Xylol or nitro dilutions which can damage the microscope.

Cleaning agents of unknown composition should be tested on an inconspicuous part of the microscope. Painted or plastic surfaces must not be tarnished or etched.

Cleaning the specimen stage

Remove bright spots on the specimen stage by wiping them off with paraffin oil or acid-free vaseline.

Cleaning glass surfaces

Remove any dust from glass surfaces with a fine, dry, and grease-free artist's hair brush, by blowing off with a bellows ball or by vacuuming.

Carefully remove obstinate dirt on glass surfaces with a clean cloth moistened with distilled water. Failing this, you may replace the distilled water with pure alcohol, chloroform or benzine.

Cleaning the objectives

Attention!

Objectives may not be disassembled for cleaning. If defects are detected on inside surfaces, send the objectives to your Leica representative for repair. We also do not advise cleaning the inside surfaces of the eyepieces.

The front lens of objectives can be cleaned as described for "Cleaning glass surfaces". The upper lens can be cleaned by blowing dust off with a bellows ball.

Removing immersion oil

Attention!

Observe the safety instructions for immersion oil!

First wipe of the immersion oil with a clean cotton cloth, then wipe over several times with ethyl alcohol.

Handling acids and alkaline solutions

Take particular care when working with acids or other aggressive chemicals.

Attention!

Always avoid direct contact of such chemicals with optical and mechanical parts.

Troubleshooting, lamp/fuse replacement

All Leica instruments are manufactured and tested with extreme care. However, if you have cause for complaint, please do not try to make any intervention on the instrument yourself. Directly contact your national agency or our central service department, the Technical Service in Wetzlar.

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Phone +49 (0) 64 41-29 28 49

Fax +49 (0) 64 41-29 22 66

Apart from preparation errors (e. g. staining or wrong specimen vessels) which are beyond the scope of this manual, there are basically two categories of defects:

Mechanical defects and
Electrical defects

Mechanical defects

Reference to possible mechanical defects has already been made in chapters "Setting up" and "Operation".

Basically, these include the improper positioning of contrast-enhancing accessories, the maladjustment of light rings or the setting of an improper condenser level.

All these possible errors have been covered in the previous chapters.

Therefore, if you are unable to obtain the desired microscope image, please read the corresponding chapters of this manual.

Electrical defects

Electrical defects may include:

1. The lamp on the microscope does not work.
2. No voltage is present.

Check the following possible causes:

The On/Off switch does not work (no illumination):

- Check that all power cords are properly connected.
- Make sure that the mains voltage is present at all power outlets used and has not been disabled with a master switch.
- After you have ruled out all possible external causes, it is possible that a fuse of the Leica DMIL microscope or the power unit is defective.

Replacing the mains fuse on the microscope

Attention!

Always disconnect the mains plug!

- Switch off the microscope.
- First disconnect the power cord from the power outlet and then disconnect it at the microscope.
- Disconnect the power cord of the power unit if required.
- Use a screwdriver to press on the latch (55.2) and remove the fuse holder (55.1).
- Remove the defective fuses from the fuse holder.
- Replace them with two new fuses of the proper type.

| Voltage rating | Type |
|----------------|------------|
| for 100 V | 2x T800 mA |
| for 115 V | 2x T800 mA |
| for 230 V | 2x T800 mA |

Note

Never use replacement fuses with a different current rating.


- Insert the fuse holder (55.1) until the latch clicks audibly in place.
- Subsequently connect the power cord (55.4) to the microscope and to the mains supply.
- If required, connect the power unit to the mains supply.

The integrated transmitted-light lamp does not work

- Make sure that the plug of the lamp cable is firmly plugged in the corresponding socket on the rear panel of the DM IL microscope stand.
- The halogen lamp may be defective.

Fig. 55

- 1 Fuse holder with voltage selector module, 2 Latch, 3 Selector module with voltage indications, 4 Mains voltage, 5 Mains connection

Replacing the 6 V/35 W halogen lamp

Attention!

Always disconnect the power plug!
Remove the protective sleeve of the lamp only after inserting the lamp. Avoid finger prints, and if there are any, always wipe them off.
Caution! Lamp and housing may be hot!

- Switch off the microscope and also the power unit, if required.
- Disconnect the power cord of the microscope and also of the power unit, if required.

- Disconnect the mains connection of the transmitted-light illumination carrier on the rear of the microscope stand (56.1).
- Use a 3 mm fixed spanner to remove the lamp housing (57.1).
- Remove the defective lamp.
- Place a new lamp (58.1) into the sockets of the lamp holder (58.2).
- Position the lamp housing and screw it in place with a 3 mm fixed spanner.
- Connect the transmitted-light illumination carrier to the mains voltage at the rear of the microscope stand (56.1).
- Connect the microscope and, if required, also the power unit to the mains supply.

The additional fluorescence lamp does not work

- Make sure that all cable connections between lamp, power unit, and mains have been properly established.
Possible causes for the failure of the fluorescence lamp can be a blown fuse of the power unit or a defective burner in the lamp housing.

Replacing the mains fuse on the power unit*

Attention!

Always disconnect the power plug first!

- Switch off the microscope and the power unit.
- Disconnect the power cords of the microscope and of the power unit.
- Remove the defective fuse from the fuse holder.


Fig. 56 Rear of microscope
1 Lamp cable connection

Fig. 57
1 Mounting screw for lamp housing

Fig. 58
1 6 V/35 W halogen lamp, 2 Lamp holder

Replacement fuses in compliance with IEC 127-2 and/or UL 198 G and/or manufacturer type:

Part no: 846-205.000-00

Designation: T 4A

Wickmann 19 195/
Schutter FST

Note

Never use replacement fuses with a different current rating.

- Connect the microscope and the power unit to the mains supply.

Replacing the 12 V/100 W halogen lamp in the 106, 107, 107/2 lamp housings

Ask a Leica field technician to show you the proper procedure for replacing the halogen lamp. All necessary steps are listed below.

Attention!

Prior to making any assembly work, always disconnect the power supply at the external transformer outlet and the microscope outlet!

- Switch off the microscope and the power unit.
- Disconnect the power cords of the microscope and of the power unit.
- Unscrew the lock screw on the microscope and remove the lamp housing.

Fig. 59 107/2 lamp housing

1 Screw for opening of lamp housing


Fig. 60 107/2 lamp housing, opened

1 Collector, 2 Lamp holder with 12 V/100 W halogen lamp



Fig. 61 106 lamp housing, opened

1 Screw for opening of lamp housing, 2 Lamp holder with 12 V/100 W halogen lamp, 3 Collector, 4 Diffuser disc

- Unscrew the screw (59.1 or 61.1) on the lid and remove the lid.
- Move the collector (61.3) to the front, if required.

Note

This step is not required for the 107/2 lamp housing.

Attention!

Leave the protective covering on the lamp until the lamp is in its holder! Avoid making finger prints or wipe them off immediately.

- Remove the defective lamp.
- Place a new 12 V/100 W halogen lamp straight into the lamp holder (60.1 or 61.2).
- Move the collector back to its original position.
- Replace the lid and secure it with the screw (59.1 or 61.1).
- Attach the lamp housing to the microscope and secure it with the lock screw.
- Connect the lamp housing to the power unit.
- Connect the microscope and the power unit to the mains supply.

Fig. 62 106 z lamp housing, opened

1 Lid, flipped up, 2 Collector, 3 12 V/100 W halogen lamp or gas discharge lamp (see Fig. 38), 4, 9 Lid mounting, 5 Reflector, 6, 8 Adjusting screws for x/y centring of reflector, 7 Focusing of reflector, 10 Mounting screws for lamp holder, 11 Socket for cut-out plug


Fig. 63 12 V/100 W lamp holder

Replacing the 12 V/100 W halogen lamp in the 106 z lamp housing*

Attention!

Prior to making any assembly work, always disconnect the power supply at the external transformer outlet and the microscope outlet!

- Switch off the microscope and the power unit.
- Disconnect the power cords of the microscope and of the power unit.
- Unscrew the lock screw on the microscope and remove the lamp housing.
- Unscrew the screws on the lid (62.4 and 62.9) with a crosstip screwdriver.
- Slightly pull out the cut-out plug from the socket (62.11) and flip up the lid.

Attention!

Leave the protective covering on the lamp until the lamp is in its holder! Avoid making finger prints or wipe them off immediately.

- Unscrew the mounting screws (62.10) on the lamp holder and remove the lamp holder (Fig. 63).
- Remove the defective lamp.
- Insert a new 12 V/100 W halogen lamp into the lamp holder.
- Insert the lamp holder and secure it with the screws (62.10).
- Plug the cut-out plug into the socket (62.11).
- Flip the lid back down and tighten the screws (62.4 and 62.9) on the lid.
- Attach the lamp housing to the microscope and secure it with the lock screw.
- Connect the lamp housing to the power unit.

Replacing the Hg and Xe lamps on the 106 z lamp housing

Attention!

- Prior to making any assembly work, always disconnect the power supply at the external transformer outlet and the microscope outlet!
- Allow the lamp housing to cool prior to opening it (at least 15 min.), explosion hazard!
- Never touch the glass parts of the burner with your hands. If required, carefully remove finger prints and dust (use alcohol if necessary).
- Adjust lamps immediately after ignition.
- Avoid frequent switching on and off because this could affect the service life and stability. Hot Hg lamps will ignite only after cooling off. It is recommended to run in new burners for a couple of hours without interruption.
- Make sure that the lamp housing is sufficiently ventilated. Never block air slots with paper etc., fire hazard!
- It is good practice to record the hours of use and to compare them to the manufacturer's specifications.
- Replace discoloured, worn burners in due time.
- We must refuse any liability for damage resulting from a possible explosion of the lamp.

- If required, disconnect the mains plug of the power unit and the microscope.
- Open the 106z lamp housing by loosening the screws (62.4), slightly pulling out the cut-out plug from the socket (62.11) and flipping up the lid of the lamp housing.
- Unscrew the safety screws (62.10) and pull out the lamp holder (Fig. 64).
- Insert the burner as follows while strictly observing the above safety instructions:
- If a plastic protective sleeve is present, leave it in place for the time being.
- Insert the burner so that the lettering is **upright** after installation. For Hg 50, Hg 100, and Xe 75, the different height of the metal base ensures installation at the proper height.
- Align any existing glass fused nipple (64.2) by rotating the burner so that the nipple is orientated **to the side** and away from the beam path.
- Insert the upper pin of the burner between the clamps of the flexible power connection and fix it with the screw (64.1) .
- Slightly unscrew the stud (64.4) in the holder.
- Insert the burner into the lower end of the metal base and retighten the stud.
- Remove the protective sleeve of the burner if it is still present.
- Place the lamp holder with the burner inserted into the lamp housing and tighten the screws (62.10).
- Close the lid of the lamp housing.
When closing the lamp housing, make sure that the pins of the cut-out plug engage in the sockets provided for this purpose.
- Retighten the screws of the lid.
- Push the cut-out plug fully in.

- Attach the lamp housing and fasten it to the microscope with the lock screw.
- Connect the lamp housing to the power unit (compare mains voltage!):

The Hg 50 W lamp is properly installed if:

1. The type is stamped on the lower socket of the lamp. The stamped lettering should be visible, i. e. **not** upside down.
2. The upper base is marked "UP".

Note:

An incorrect installation will reduce the lamp brightness to about 60% and will considerably limit the useful life of the lamp.

• **Attention:**

Make sure that the markings on the lamp base and on the power unit are the same. For example, if the lamp base is marked L1, L1 must also be set on the power unit to make full use of the lamp and not to shorten its life.



• **Important:**

Make sure to dispose of worn burners in an environmental-friendly manner.

Fig. 64 Lamp holders for gas discharge lamps

1 Upper clamp, 2 Fused nipple of burner, 3 Lower clamp, 4, 6 Mounting holes for holder, 5 Sockets for cut-out plug, 7 Protective sleeve

Storage

Protect your microscope from dust after use by covering it with the protective cover.

Store the microscope in a cabinet in which the temperature is ≥ 5 °C above room temperature. The cabinet must be equipped with ventilation openings which are filled e. g. with cotton swabs to protect the interior from dust. If such a storage is not possible, place the microscope into a closed container which contains a desiccant (e. g. silica gel) .

Packaging and transport

Ship or transport the microscope and its accessories in its original packing. A delivery note with all necessary information must be included in the shipping container.

Technical description

Due to basic physical principles and the physiology of the human eye, all imaging techniques, not only the microscope, are subject to limitations in performance. For proper use of the microscope you should therefore know and observe the following information.

Performance data of objectives

The Leica DM IL microscope is based on a tube length of ∞ (infinite) and a focal length of the tube lens of $f = 200$ mm.

Attention:

For this reason, only objectives with the engraving ∞ and M 25 thread may be used.

Objective labelling

Examples and meaning of the symbols:

$\infty / -$
C PLAN 10x/0.22

$\infty / 0.17$
C PLAN 40x/0.65

$\infty / 0 / D$
N PLAN 50x/0.75

∞

Objective for infinite tube length (∞).

—

The objective can be used **with and without** a coverglass.

0.17

The objective may only be used with a coverglass of the standard 0.17 mm thickness. Use without a coverglass or with a coverglass of a very different thickness will result in a distinct drop in performance, especially for objectives with high apertures (see below).

0

Use **without** coverglass, e.g. for cell smear specimens, incident light. Not suitable for inverse microscopes.

D (or A, B, C)

Pupil position of the objective (important e.g. for Integrated Modulation Contrast IMC).

Objective type (performance class):

| | |
|--------------------|--|
| C Plan | Semi-Planachromat |
| N Plan | Planachromat |
| PL FLUOTAR® | Plan-Semiapochromat |
| PL APO | Planapochromat |
| HC | Harmonic Components |
| X | Universal application potential, also backwards-compatible with Delta optics (= predecessor of the HC optics). |
| L | Long working distance. |
| 10 x/0.22 | Magnification and aperture. The aperture (pick-up angle) determines resolution, field depth, contrast and brightness. Objectives with a built-in iris diaphragm are engraved with their maximum and minimum aperture, e. g. 0.85 – 0.55. ! Attention: Objectives with a built-in iris diaphragm! The knurled ring may only be used to adjust the diaphragm, not for screwing in and out. Risk of damage! |
| OIL, W, IMM | Immersion objectives for: oil, water, universal (oil, water, glyzerine, etc.) |
| PH | PH = Phase contrast objective, the corresponding light ring in the condenser is also indicated, e. g. PH2. |

| | |
|---------------|---|
| BD | BD = Brightfield/Darkfield; objectives for incident light microscopy with M 32 objective thread. |
| P, POL | Strain-free objective for quantitative polarised light microscopy. |
| U-V-I | Special apochromatic correction, i.e. parfocal from U ltraviolet through V isual to near I nfrared (from about 340 nm to 1000 nm). |

Colour codes of objectives

In compliance with DIN/ISO standards, the magnification of each objective is indicated by a colour ring:

| | | | | | | | | | |
|------------------------------|-----------|------------|------------|-------------|--------|--------|----------|-------|------|
| 100x 125x 150x 160x | 63x | 40x 50x | 25x 32x | 16x 20x | 10x | 6.3x | 4x 5x | 2.5x | 1.6x |
| white | dark blue | light blue | dark green | light green | yellow | orange | red | brown | grey |

Immersion objectives are also marked with a second, lower colour ring:

| | |
|---------------|--|
| black | Oil or Imm (= Universal objective for oil, water, glyzerine) |
| white | Water |
| orange | Glyzerine |

Performance data of eyepieces

Our product range includes the following eyepieces for the Leica DM IL:

Eyepieces for viewing with DM ILB or DM ILT tubes

| Magnification/ Field of view | Eyepiece port +) | |
|---------------------------------|---------------------|---|
| 10 x/18 | | |
| 10 x/18 | | M |
| 10 x/20 | | |
| 10 x/20 | | M |
| 15 x/14 | | |

Eyepiece tube diameter: 23.2 mm

Eyepieces for viewing with tubes from the DM L range

| Leica eyepiece type | Magnification/ Field of view | Eyepiece port +) | |
|------------------------|---------------------------------|---------------------|----|
| HC PLAN | 10 x/20 | | M |
| HC PLAN | 10 x/20 | | |
| HC PLAN | 12.5 x/16 | | M |
| HC PLAN | 10 x/20 | | MF |

Eyepiece tube diameter: 30 mm

- +)
- = with detachable or fold-back glare protection, for use with or without glasses
 - M = adjustable eyelens (dioptre compensation) and mount for graticules of 19 mm or 26 mm for HC eyepieces
 - MF = with illuminated graticule

The eyepiece type LEITZ PERIPLAN® must not be used!

Eyepieces of the earlier type L PLAN may only be used with earlier-type tubes (before about 1998) **without** engraving HC!

Eyepiece field of view number

For a certain microscope configuration a certain eyepiece field of view number must not be exceeded (see below), e. g. 20. If the maximum field of view is exceeded, this may result in a disturbing loss of definition or vignetting at the edge of the image, → following pages!

The eyepiece field of view (fov) designates the diameter of the intermediate image in the eyepiece in mm, i. e. the diameter of the circular diaphragm which limits the image format and which lies within the eyepiece.

This fov is indicated on the eyepiece after the magnification, e. g. 10x/20.

A maximum fov **of 20** is recommended for the Leica DM IL microscope.

The maximum admissible eyepiece field of view number of a certain configuration is derived from the following instrument data:

| | |
|--|------------------|
| Field performance of objectives | see below |
| Field performance of intermediate module(s) | see below |
| Field of view number of tubes | → p. 74 |
| Condenser characteristics | → p. 75 |

The decisive value is always the **smallest**.

If e. g. the intermediate modules (see below) only allow the field of view number 20 while the objectives and tube allow 25, the maximum admissible fov number for the eyepiece is 20. Eyepieces with a fov number of 25 can cause vignetting. In detail, the following applies:

The diameter of the **observed object field** is calculated by dividing the field of view diameter by the objective magnification and the magnification of the microscope optics.

Example:

Eyepiece 10x/20

Objective PLAN 4/0.10

Magnification factor of Leica DM IL microscope optics 1x

Observed object field

$$\frac{20 \text{ mm}}{4 \times 1} = \varnothing 5 \text{ mm}$$

The **overall magnification** of the microscope is calculated by multiplying the magnification of the eyepiece with the magnification of the objective and the magnification of the microscope optics.

Example:

Eyepiece 10x/ 20

Objective PLAN 4/0.10

Magnification factor 1x

Overall magnification $10 \times 4 \times 1 = 40x$

Field performance of objectives

The engraving on the objectives does not indicate their field performance. This performance can vary slightly within the same class of objectives, e.g. the lower objective magnifications may well have slightly higher values than the approximative values given below:

| Objective series | max. recommended eyepiece fov | | | |
|------------------------|-------------------------------|------------|----|----|
| | 15 | 20 | 22 | 25 |
| Achromats | ██████████ | | | |
| C PLAN Achromats | ██████████ | | | |
| APO L Apochromats | ██████████ | | | |
| N PLAN Planachromats | ██████████ | ██ | | |
| PL FLUOTAR® Semiapo. | ██████████ | ██████████ | | |
| PL APO Planapochromats | ██████████ | ██████████ | ██ | |

A current data sheet covering all Leica objectives can be obtained from your local Leica representative.

Field performance of intermediate modules

The maximum admissible field performance of the intermediate modules is derived from the type designation listed in the following table and also on your invoice. Each type designation consists of two values which are separated by a slash, e.g. Ergo module **L 2/25**.

The first value (2 in the example) is a relative measure (height index) of the overall height of the module. Multiplying the height index by the factor 15 yields the amount in **mm** by which the viewing port or the overall height of the microscope is increased, i.e. $2 \times 15 = 30 \text{ mm}$. The second value (25 in the example) is the maximum field of view number which is possible with this module.

Ergo module L 2/25

Magnification changer L 3/25

Drawing facility L 3/20

Discussion facility L 3/20 (2 viewers)

Performance data of filters

| Filter | Application |
|--------------------------------------|---|
| Grey filter N/Neutral filter | Grey filters (neutral filters) are used to attenuate the light without influencing the colour temperature. The engraved value, e.g. N16, indicates the attenuation. N16 indicates a reduction to $1/16 = 100/16 = 6.25\%$ transmission. |
| Green filter, GR panchromatic | Contrast enhancement for b/w images. |
| DLF | Conversion filter (daylight filter blue, similar to CB12) for colour photography with daylight film, integrated in filter magazine. |
| BG38 (Blue filter) | Suppression of red in fluorescence applications (integrated in fluorescence illuminator). |
| ALF | Artificial light filter for colour photography with artificial light film, to enhance colour contrast. |
| BG20 | For red enhancement in Polaroid images. |
| VG9 (Green filter) | Contrast enhancement for chromosome photography. |
| CB1.5, CB3 | Conversion filter blue: To increase the colour temperature when using special lamps. |
| CR1.5 | Conversion filter red: To reduce the colour temperature, e.g. from 6000 K (colour temperature of a Xe lamp) to 5500 K (colour temperature for daylight film). |
| BG23 | Contrast enhancement of the complementary colours blue and red on b/w film. |

Performance data of tubes

The tube changer is the same as with the upright microscope stands.

The tubes can be rotated and changed.

DM ILB binocular tube

The binocular tube consists of a body with the tube changer ring attached to the lower side. The tube lens has a factor of 1x. The Siedentopf binocular tube permits to adjust the interpupillary distance from 55 mm to 75 mm while maintaining the tube length. The viewing angle is 45°. The tube has an adjustable eyepiece tube. It offers a field of view number of 20. When intermediate tubes are used, the maximum possible field of view number is 18 (e. g. when using the drawing facility).

DM ILT trinocular tube

The trinocular tube consists of a body with the tube changer ring attached to the lower side. The tube lens has a factor of 1x. The Siedentopf binocular tube permits to adjust the interpupillary distance from 55 mm to 75 mm while maintaining the tube length. The viewing angle is 45°. The tube has an adjustable eyepiece tube. It offers a field of view number of 20. When intermediate tubes are used, the maximum possible field of view number is 18 (e. g. when using the drawing facility or the multi-discussion facility).

The lateral documentation port is to be fitted with HC components only.

The tube contains a switchable mirror with two positions:

- a) 100 % of light directed to binocular tube
- b) 100 % of light directed to photo tube

The optical axis of the documentation port is offset to the left by 88 mm. This offers a free view of the specimen when using this tube.

Fig. 65 DM ILB binocular tube

Fig. 66 DM ILT trinocular tube

Tubes from the DM L range

The use of tubes from the DM L range requires a tube adapter. The DM IL/L tube adapter is an intermediate tube with a length of 60 mm and without optics which is used to adapt the pupil position. It has a tube changer ring at the bottom and a tube changer area for the DM L tubes at the upper end.

The following tubes can be used:

Binocular tube HC LB 0/3/4 and HC LBP 0/3/4

Trinocular tube HC L1T 4/5/7 and HC L1TP 4/5/7

Trinocular tube with 3 switch positions HC L3TP 4/5/7

Ergo tube, binocular HC LVB 0/4/4

Ergo photo tube, trinocular HC L1VT 0/4/4

Field of view number of DM L tubes

The type designation of DM L tubes also contains a combination of numbers to indicate the maximum admissible eyepiece fov number, e. g. binocular tube HC LB **0/3/4**.

The numbers **0/3/4** indicate the maximum admissible height index of the intermediate modules for the eyepiece field numbers **25, 22 and 20**.

In the above example, this means:

1st number (**0**): fov 25 can only be obtained if the tube is directly mounted onto the microscope, i. e. **without** an intermediate system.

2nd number (**3**): fov 22 is only possible up to height index 3, e. g. magnification changer **L3/25**.

3rd number (**4**): fov 20 up to max. height index **4**, e. g. **2** Ergo modules **L 2/25**.

If the number is replaced with a hyphen –, e. g. monocular tube LMP –/–/7, the tube cannot be used for the corresponding fov at all; in the above example, this means it cannot be used for fov 25 and 22, while fov 20 is possible up to index 7.

Exceeding the admissible values can result in vignetting (shading at the edges of the image).

HC label: Only eyepieces of type **HC** PLAN and widefield 16x and 25x can be used. If there is no HC label, this means that eyepieces of the type Leica **L** PLAN have to be used.

Additional examples:

0/4/4 Field of view 25 is only possible with direct tube mounting to the microscope stand (height index of intermediate modules = 0), provided that suitable objectives are used.

Fov 20 and 22 are possible up to height index 4, e. g. with the fluorescence device. The addition of a further module would not be admissible; a solution to the problem would be a tube with the following characteristics:

4/5/7 Fov 25 is possible up to height index 4 (e.g. 2 Ergo modules L2/25 or magnification changer L3/25). Fov 22 is possible up to height index 5, fov 20 up to height index 7 (e. g. illuminator LRF 4/20 plus magnification changer L3/245).

–/–/7 The tube allows fields of view up to 20 mm only. If intermediate modules are used, the total of their height values must not exceed 7.

Performance data of condensers

S 90 0.23 condenser

Suitable for lab vessels up to a height of 90 mm and objectives with numerical apertures up to 0.50. Offers an optimal matching of light and phase rings up to liquid levels of 35 mm in phase contrast applications.

S 55 0.35 condenser

Suitable for lab vessels up to a height of 55 mm and objectives with numerical apertures up to 0.60. Offers an optimal matching of light and phase rings up to liquid levels of 50 mm in phase contrast applications.

Applications of condensers S 90 and S 55:

| Illumination method | S 90 Objectives | Light rings/ Accessories | S 55 Objectives | Light rings/ Accessories |
|--------------------------------|------------------|--------------------------|---|-------------------------------|
| Bright-field | 25x 4x – 100x | (with diffuser) – | 2.5x 4x – 100x | (with diffuser) – |
| Phase contrast | 5x 10x – 20x | Phaco 0 Phaco 1 | 5x 10x – 20x 40x – 63x | Phaco 0 Phaco 1 Phaco 2 |
| Integrated Modulation Contrast | not possible | | C PLAN 10x/0.22 AP 32.2 C PLAN L 20x/0.30 D C PLAN L 40x/0.50 D and all objectives with pupil position D | IMC slide |

Incident-light fluorescence illumination*

To obtain a brighter image for incident-light fluorescence applications, the Leica DM IL microscope is preferably fitted with mercury and Xenon gas discharge lamps but can also be operated with a 12 V/100 W halogen lamp.

Performance data of lamp housings

106 lamp housing*

The 106 lamp housing is fitted with a 12 V/100 W halogen lamp. The lamp holder can be centred in the x/y plane. The aspherical collector can be focused. The 106 lamp housing does not have a reflector but is fitted with a diffuser disc and a heat-absorbing filter.

106 z lamp housing*

Same as 106 lamp housing but with centrable and focusable reflector and 4- to 6-lens collector. A quartz collector is available on request. The following lamps can be used, each with the corresponding special lamp holder:

- Halogen lamp 12 V/100 W
- Hg ultra-high pressure lamp 50 W (AC)
- Hg ultra-high pressure lamp 100 W (DC, stabilised)
- Hg ultra-high pressure lamp 100 W (DC, stabilised Type 103 W/2)
- Xe ultra-high pressure lamp 75 W (DC, stabilised)

107/2 lamp housing

The shield connection of the lamp housing is screwed down to the equipotential bonding point on the 12 V/100 W power unit. This lamp housing for incident-light and transmitted-light operation has a 1-lens fixed collector and a fixed 12 V/100 W lamp.

Note:

The LH 105 lamp housings have been replaced with the LH 106 lamp housings. However, they are compatible to the LH 106 lamp housings and can also be used.

| Type | Typical lifetime |
|--|------------------|
| Halogen lamp 6 V/35 W | 50 h |
| Halogen lamp 12 V/100 W | 50 h |
| Hg ultra-high pressure lamp 50 W (AC) | 100 h |
| Xe high pressure lamp 75 W (DC) | 400 h |
| Hg ultra-high pressure lamp 100 W (DC) | 200 h |
| Hg ultra-high pressure lamp 100 W (DC, Type 103 W/2) | 300 h |

Overview of lamp housings with order numbers

Non-centrable lamp housings

| | LH 106 | LH 107, left | LH 107/2 | LH 35/2 |
|---|--------------------|--------------|------------------------|---------|
| 6 V/35 W 12 V/100 W, 0.55 m 12 V/100 W, 2.0 m 12 V/100 W, 2.0 m, shielded | 504 058 504 059 | 504 086 | 504 080 504 085 | 504 088 |

Centrable lamp housings

| | L 106, right | | LH 106, left |
|--|--|-----------------------------------|------------------------|
| | 4-lens | 6-lens | 6-lens |
| 12 V/100 W, 0.55 m 12 V/100 W, 2 m 12 V/100 W, 2.9 m | 507 070 504 071 | 504 087 504 065 | |
| Hg 100 W, w. cable Hg 100 W, w. cable, 3 m Hg 100 W, w/o cable Hg 50 W Xe 75 W | 504 068 504 069 504 083 504 066 | 504 062 504 063 504 061 | 504 090 504 089 |

General specifications

For indoor use only

| | |
|------------------------|----------------------|
| Supply voltage: | 100/115/230 V~ |
| Mains frequency: | 50 – 60 Hz~ |
| Power consumption: | 50 VA |
| Fuses: | 2x T800 mA |
| Operating temperature: | 10 – 36 °C |
| Relative humidity: | 0 – 80 % up to 30 °C |
| Overvoltage category: | II |
| Pollution degree: | 2 |

Specifications of the power unit

General specifications

For indoor use only

| | |
|-----------------------|----------------------|
| Supply voltage: | 90 – 250 V~ |
| Mains frequency: | 50 – 60 Hz |
| Power consumption: | 160 W |
| Fuses: | T 4 A |
| Ambient temperature: | 10 – 36 °C |
| Relative humidity: | 0 – 80 % up to 30 °C |
| Overvoltage category: | II |
| Pollution degree: | 2 |

Specifications

| | |
|---|---|
| Lamp DC voltage: | Adjustable from 2.5 V ± 5 % to 12 V – 5 %/8.5 A |
| Voltage setting: | Potentiometer 5 Kohms Turn clockwise for maximum intensity |
| Maximum lamp voltage: | 12.0 V for 90 V to 250 V~ |
| Soft start: | Rise time to maximum output voltage 0.2 to 1 second |
| Mains voltage dependence | |
| U_N = Mains voltage | |
| U_{La} = Lamp voltage | |
| U_N : 90 – 250 Vac, U_{La} = 12 V: | < – 5 % |
| U_N : 90 – 250 Vac, U_{La} ≤ 11 V: | < ± 1 % |
| U_N : 100 – 130 Vac, U_{La} ≤ 11 V: | < ± 0.5 % |
| U_N : 200 – 250 Vac, U_{La} ≤ 11 V: | < ± 0.5 % |
| Lamp voltage drift 0 to 10 min: | < 2 % |
| Efficiency: | approx. 75 % |
| Short-circuit and open-circuit proof | |
| Life expectancy: | > 5.0000 h |

Main wear and replacement parts

| Order No. Part No. | Component | Used for |
|---|--|--|
| <u>Spare lamps</u> | | |
| 500 322 | Halogen lamp 6 V/35 W OSRAM 64275 | Integrated illumination transmitted light |
| 500 974 | Halogen lamp 12 V/100 W | Lamp housing 105 |
| 500 137 | Hg ultra-high pressure lamp 50 W | Lamp housing 106 z |
| 500 138 | Hg ultra-high pressure lamp 100 W | Lamp housing 106 z |
| in preparation | Hg ultra-high pressure lamp 100 W (103 W/2) | Lamp housing 106 z |
| 500 139 | Xenon high pressure lamp 75 W | Lamp housing 106 z |
| <u>Tools/Adjustment keys</u> | | |
| 016-500.020-001 | Hexagonal screwdriver | Assembly and adjustment |
| 020-434-045 | 2.5 mm Allen key, angled, shortened | Assembly of heating stage and illumination mirror |
| <u>Screw covers for unused nosepiece positions</u> | | |
| 020-422.570-000 | Screw cover M25 | Objective nosepiece |
| <u>Spare eyecups (glare protection) for HC PLAN eyepieces</u> | | |
| 021-500.017-005 | Eyecup for HC PLAN | Eyepiece 10x/25 |
| 021-264.520-018 | Eyecup for HC PLAN | Eyepiece 10x/22 |
| 021-264.520-018 | Eyecup for HC PLAN | Eyepiece 10x/20 |
| 021-252.505-012 | Eyecup for standard eyepiece | Eyepiece 10x/18 |
| 004-168.001 and 120 | Eyecup for standard eyepiece | Eyepiece 10x/20 |
| 004-168.001 and 120 | Eyecup for standard eyepiece | Eyepiece 10x/18 M |
| <u>Immersion oil</u> in compliance with DIN/ISO, fluorescence-free | | |
| 513 787 | 10 ml | Objectives OIL and IMM |
| 513 522 | 100 ml | |
| 513 788 | 500 ml | |
| <u>Spare fuses</u> in compliance with IEC 127-2 and/or UL 198 G and/or manufacturer type: | | |
| 846-205.000-00 | T 4A Wickmann 19 195/ Schutter FST | Leica 12 V/100 W power unit |
| 826-252.000-00 | T 800 mA Wickmann 19 195/ Schutter FST | Leica DM IL microscope |
| <u>Ignition capacitor</u> | | |
| 302-053.023-001 | Ignition capacitor | Power unit HG 50 (500 277) |

EU Declaration of conformity

We hereby declare that the device described below, both in its basic design and construction and in the version marked by us, conforms to the relevant safety- and health-related requirements of the appropriate EU directives.

This declaration shall cease to be valid if modifications are made to the device without our approval.

Product: **DM IL/DM ILM**

Model: **Microscope**

Identification no.: 301-135.001
 301-135.002
 301-135.003

EU directives: Low voltage: 73/23/EWG
 EMC: 89/336/EWG


Harmonised standards applied: EN 61 010-1/1993
 EN 50 081-1/1992
 EN 50 082-1/1997

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Copyright © Leica Microsystems Wetzlar GmbH · Ernst-Leitz-Straße · 35578 Wetzlar · Tel. (064 41) 29-0 · Fax (064 41) 29-25 99 LEICA and the Leica logos are registered trademarks of Leica Technology BV.
Leica Microsystems Wetzlar GmbH Tel. +49 (0) 64 41-29 0